




BEACH SET


SIZE:



BASIC DIMENSIONS:



LENGTH WITH 12" FRONT WHEEL = 143 cm




FRAME WIDTH = 63 cm

WHEEL TRACK MAX = 110 cm


DIMENSIONS FOR TRANSPORT:

*) LOWER PART OF BODY MAX LENGTH

***) UPPER PART OF BODY MAX LENGTH




LENGTH = 118 cm



WIDTH MAX = 63 cm

A - A



36 cm




BEACH SET

SIZE:




BASIC DIMENSIONS:



LENGTH WITH 12" FRONT WHEEL = 167 cm

HEIGHT = 95 cm
HANDLE HEIGHT MAX - STANDARD = 110 cm
LONG = 117 cm




FRAME WIDTH = 63 cm
WHEEL TRACK MAX = 110 cm

16 cm


DIMENSIONS FOR TRANSPORT:

- *) LOWER PART OF BODY MAX LENGTH
- **) UPPER PART OF BODY MAX LENGTH




LENGTH = 138 cm

HEIGHT = 64 cm



WIDTH MAX = 63 cm

A - A



36 cm


ROVER

BEACH SET


SIZE:



BASIC DIMENSIONS:



HEIGHT = 103 cm
HANDLE HEIGHT MAX - STANDARD = 112 cm
LONG = 119 cm





FRAME WIDTH = 63 cm

WHEEL TRACK MAX = 110 cm

LENGTH WITH 12" FRONT WHEEL = 191 cm


DIMENSIONS FOR TRANSPORT:

- *) LOWER PART OF BODY MAX LENGTH
- **) UPPER PART OF BODY MAX LENGTH



WIDTH MAX = 63 cm

A - A



36 cm