




BASIC SET


SIZE:



BASIC DIMENSIONS:




LENGTH WITH 6" FRONT WHEEL = 118 cm




DIMENSIONS FOR TRANSPORT:

*) LOWER PART OF BODY MAX LENGTH

**) UPPER PART OF BODY MAX LENGTH




LENGTH = 118 cm



WIDTH MAX = 63 cm

A - A



36 cm




BASIC SET

SIZE:




BASIC DIMENSIONS:



LENGTH WITH 6" FRONT WHEEL = 138 cm

HEIGHT = 95 cm
STANDARD = 110 cm
HANDLE HEIGHT MAX - LONG = 117 cm




FRAME WIDTH = 63 cm
WHEEL TRACK MAX = 76 cm

18 cm


DIMENSIONS FOR TRANSPORT:

- *) LOWER PART OF BODY MAX LENGTH
- **) UPPER PART OF BODY MAX LENGTH




LENGTH = 138 cm

HEIGHT = 64 cm



WIDTH MAX = 63 cm

A - A



36 cm

33 cm



42 cm



BASIC SET




SIZE: **L**

BASIC DIMENSIONS:



DIMENSIONS FOR TRANSPORT:

- *) LOWER PART OF BODY MAX LENGTH
- **) UPPER PART OF BODY MAX LENGTH






JOGGING SET

SIZE:




BASIC DIMENSIONS:



LENGTH WITH 14" FRONT WHEEL = 149 cm

HEIGHT = 87 cm
HANDLE HEIGHT MAX - STANDARD = 108 cm
LONG = 115 cm




FRAME WIDTH = 63 cm
WHEEL TRACK MAX = 76 cm

18 cm


DIMENSIONS FOR TRANSPORT:

- *) LOWER PART OF BODY MAX LENGTH
- **) UPPER PART OF BODY MAX LENGTH




LENGTH = 118 cm

HEIGHT = 72 cm



WIDTH MAX = 63 cm

A - A



36 cm

32 cm

42 cm


ROYER

JOGGING SET

SIZE:


M

BASIC DIMENSIONS:



LENGTH WITH 14" FRONT WHEEL = 172 cm

HEIGHT = 95 cm
STANDARD = 110 cm
HANDLE HEIGHT MAX - LONG = 117 cm




FRAME WIDTH = 63 cm
WHEEL TRACK MAX = 76 cm

18 cm


DIMENSIONS FOR TRANSPORT:

- *) LOWER PART OF BODY MAX LENGTH
- **) UPPER PART OF BODY MAX LENGTH




LENGTH = 138 cm

HEIGHT = 64 cm



WIDTH MAX = 63 cm

A - A



33 cm

42 cm



36 cm

ROVER

JOGGING SET

SIZE: **L**

BASIC DIMENSIONS:



DIMENSIONS FOR TRANSPORT:

- *) LOWER PART OF BODY MAX LENGTH
- **) UPPER PART OF BODY MAX LENGTH

